

HHD Inventory Manager

User Manual - Version 1.1.0

HHD Inventory Manager is a Windows application for tracking home hemodialysis supplies, treatments, received inventory, alerts, and historical inventory levels.

Dashboard and Responsive Tables

Inventory tables automatically adjust row height and column width to the selected font size and to the space available inside each panel. The item-name and notes columns receive the most flexible space so the surrounding panels remain intact.

- Use the minus and plus controls above inventory tables to adjust the inventory font size.
- Green Alert means OK, flashing yellow means LOW, and flashing red means RE-ORDER.
- NxStage and DaVita dashboard panels remain evenly divided.

Themes

Open Settings / Items and choose an Application Theme. Available themes are Medical Blue, Beige, and Dark. The selected theme is saved and restored when the application opens.

Theme	Description
Medical Blue	Original medical-blue interface.
Beige	Warm light interface with dark text and high-contrast controls.
Dark	Neutral dark interface with light text.

Treatments and Treatment Calendar

The main left-panel button is labeled Treatments. Use it to record regular, extra, missed, or incomplete treatments and to review or edit treatment notes.

Treatment Calendar provides Month and Week views. Use Previous, Today, and Next to navigate.

- Green day: treatment performed.
- Yellow day: incomplete treatment.
- Red day: missed treatment.
- Unmarked days retain the normal calendar background.

Database Import and Export

Import Database restores a selected SQLite database after creating a safety backup. Export Database creates a portable copy of the current database at a location selected by the user.

Inventory History

Inventory History charts calculated inventory over 1 week, 1 month, 3 months, 6 months, 1 year, or all available time. The horizontal axis always spans the selected period. Dates before an item's baseline are marked as unavailable rather than estimated.

SAK Usage Rule

For SAK, regular, extra, and incomplete treatments each consume one reusable-session slot. Missed treatments consume no SAK slot. With a SAK configured for two sessions, each attempted treatment deducts 0.5 SAK.

Data and Backups

The live database and settings are stored in the current user's AppData folder. Rolling backups are stored in Documents/HHD Inventory Backups.

Created by Eduardo A. de Carvalho, husband and caregiver of Joelle.