

HHD Inventory Manager

User Manual - Version 1.0.0

HHD Inventory Manager is a Windows desktop application for tracking Home Hemodialysis supplies. It is designed to help maintain a clear inventory record for NxStage and DaVita supplies, calculate estimated usage, and highlight supplies that need reorder attention.

Important safety note: This software is an inventory tracking tool only. It does not replace medical, DaVita, NxStage, physician, nurse, or clinical instructions. Always verify physical inventory before treatment.

Contents

1. Overview
2. Installation and first launch
3. Dashboard
4. Settings and schedule
5. Item configuration
6. Received inventory
7. Dialysis session log
8. Current inventory calculations
9. Backups, import, and moving to another computer
10. CSV export
11. Recommended workflow
12. Troubleshooting

1. Overview

The program organizes supplies into two inventory groups: NxStage Supplies and DaVita Supplies. Each item has configurable usage rules, thresholds, and baseline inventory information. The application stores data locally in a SQLite database.

Inventory Group	Typical Items
NxStage Supplies	SAK, PAK, Cartridge, Dialysate Hanging Bags, Warmer Lines
DaVita Supplies	Heparin, 10CC Syringe, Medical Gloves, Syringe Needles, 15 Gauge Tulip Needles, Alcohol Pads, Iodine Pads, Saline bags, Mr. Peanut connector, gauze, Paper Towels, Chloramine Test Strips

2. Installation and First Launch

If using the installer, run the setup program and follow the prompts. If running from source, install Python 3.11 or newer and run `run_HHD_Inventory_Manager.bat`.

On first launch, the program creates a local database named `hhd_inventory.db`. It also creates a settings file when window settings or preferences are saved.

3. Dashboard

The Dashboard is the main overview screen. It displays the patient name, local date and time, inventory panels for NxStage and DaVita supplies, and at-a-glance status cards.

- Units Left shows the current calculated inventory.
- Weeks Left estimates how long the inventory may last based on item settings.
- Status shows OK, LOW, or RE-ORDER.
- Re-order and Low counts show how many items need attention.

4. Settings and Schedule

Open Settings / Items to configure the patient name, dialysis sessions per week, and the week's first session day. The sessions per week value is used when forecasting weeks remaining.

Example: if sessions per week is set to 4 and an item uses 5 units per session, the forecasted session-based weekly usage is 20 units per week.

5. Item Configuration

Each inventory item can be edited. Double-click an item or select it and choose Edit / Rename Selected. You can also add new items to either inventory group or remove items from active inventory.

Field	Meaning
Item Name	The name shown in the inventory lists.
Inventory Group	NxStage Supplies or DaVita Supplies.

Field	Meaning
Current/Baseline Units	The inventory count on the baseline date.
Inventory Count Date	The date the baseline count was taken.
Units Used Per Session	How many units the item normally consumes per dialysis session.
Additional Units Used Per Week	Weekly usage not tied to session count, such as one medical tape per week.
Reusable For # Sessions	How many sessions one unit lasts. Example: SAK may be reusable for 2 sessions.
Lifespan Days Per Unit	For items like PAK that last a number of days instead of a per-session rate.
Low Threshold	When units reach this value, status becomes LOW.
Re-order Threshold	When units reach this value, status becomes RE-ORDER.
Auto-calculate usage	When enabled, logged sessions reduce calculated inventory for that item.

6. Received Inventory

Use Add Received Inventory when new supplies arrive. Select the item, enter the received date, units received, and optional notes. The program stays on the Received Inventory screen so multiple items can be entered without returning to the dashboard.

- Use units, not boxes.
- Enter the actual received date if different from today.
- Received inventory is added to the current calculated inventory.
- Recent received entries are shown on the same screen.

7. Dialysis Session Log

Use Log Session to record dialysis sessions. Session entries are used to calculate inventory usage for items that have automatic session usage enabled.

Session Type	Suggested Equivalent
Regular Session	1
Extra Session	1
Missed Session	0
Incomplete Session	0.5 or another appropriate decimal value

8. Current Inventory Calculations

The program calculates current inventory from the baseline count, received inventory, corrections, and calculated usage.

Current units = baseline units + received units + corrections - calculated usage

Calculated usage can come from logged sessions, weekly usage, reuse rules, and lifespan rules. For forecasting weeks remaining, the configured sessions per week are used.

If an item should be tracked manually or used only in emergencies, disable automatic session usage or set usage fields accordingly.

9. Backups, Import, and Moving to Another Computer

The application automatically backs up the database every 10 minutes. The backup is a rolling backup that overwrites the previous backup to avoid creating many files.

File	Purpose
hhd_inventory.db	The live working database.
Documents\HHD Inventory Backups\HHD_Inventory_Backup_Current.db	Rolling automatic backup of the database.

File	Purpose
Documents\HHD Inventory Backups\HHD_Settings_Backup_Current.json	Rolling automatic backup of local settings.
Documents\HHD Inventory Backups\HHD_Inventory_Backup_Before_Import.db	Safety backup made before importing another database.

To move to another computer, copy the database file or the rolling backup file to the new computer and use Import Database from the program.

10. CSV Export

Use Export CSV to create a spreadsheet-compatible inventory summary. The export includes item names, current units, status, weeks remaining, baseline details, received totals, and calculated usage.

11. Recommended Workflow

- Set patient name and dialysis schedule first.
- Review each item and configure units per session, reuse count, weekly usage, and thresholds.
- Enter initial baseline inventory counts using units only.
- Use Add Received Inventory whenever supplies arrive.
- Use Log Session after each dialysis session or when an extra, missed, or incomplete session occurs.
- Review the Dashboard for LOW or RE-ORDER status.
- Verify physical inventory regularly and update baseline or item settings when necessary.

12. Troubleshooting

Problem	Suggestion
Inventory looks wrong	Check baseline count, baseline date, logged sessions, received inventory, and item usage settings.
An item should not auto-deduct	Edit the item and disable automatic session usage or set usage values to zero.
Program moved to another computer	Use Import Database to load the prior hhd_inventory.db or backup database.
Taskbar icon does not update	Rebuild the EXE using build_exe.bat and launch the built EXE rather than running the Python script directly.
Need a backup	Copy hhd_inventory.db or the rolling backup file from Documents\HHD Inventory Backups.